

Acne Treatment Results



Photos courtesy of J. Matthew Knight, M.D.



Photos courtesy of Chehao Hsu, M.D.



Photos courtesy of Chehao Hsu, M.D.




Photos courtesy of Yu-Pin Peng, M.D.

PD Dr. med. Marianne Placzek-Neukirchen
Statenkliniek
Dermatologie
Frankenslag 357
NL-2582 HP Den Haag
Tel. +31 70 3220 920
Fax. +31 70 3220 921
info@statenkliniek.nl
www.statenkliniek.nl



WWW.AESTHETIPEDIA.COM

©2017 All Rights Reserved. The Lumenis Group of Companies



Acne Solution with M22™

Transform your skin



Why treat acne with M22™?

It's a fast "lunch-time" treatment with minimal downtime

A great alternative for oral medication

It improves your skin appearance and helps you regain your self-confidence after just a few treatments

What is it?

The Lumenis® acne solution

The Lumenis® acne solution enables your physician to greatly improve your skin appearance and reduce your acne, while also treating skin redness and pigmentation.

The unique Lumenis IPL technology on the M22 system applies light flashes over the treated area, triggering a response that will eventually eliminate the bacteria within the pores. The M22 system offers a non-invasive solution for inflammatory acne and offers a great alternative for oral medication.

Is it right for me?

Lumenis energy-based treatments are not suitable for everyone and carry some risks. IPL for acne is not suitable for those who suffer from active viral or fungal skin conditions, or skin cancer. Risks may include: redness, swelling and change of pigmentation. Be sure to consult with your treatment provider before choosing this treatment.

What should I expect?

When will I see results?

During the Lumenis M22 acne treatment, short pulses of intense light are emitted from an applicator onto your skin. The light penetrates through the most superficial skin layers and helps eliminate acne bacteria and reduce the inflammation or excess sebum (skin oil) production that characterizes acne. Depending on acne severity and skin type, it may typically take 4-6 sessions with an average recess of 2-4 weeks in between treatments to see improvement.

What can I expect after treatment?

Immediately following treatment, you may experience some redness, depending on your customized treatment settings. Daily activities can be resumed the very same day.




How does it work?

M22 IPL acne technology

IPL is one of the most commonly administered form of light therapy in dermatological applications. IPL systems are designed to deliver many rapid, highly intense and controlled pulses of light – the controlled pulsing prevents thermal damage to your skin and minimizes discomfort. The M22 technology treats your acne with specific wavelengths of light targeting bacteria in the skin, as well as inflamed sebaceous glands that contribute to break-outs. IPL energy (at various wavelengths) is also used to selectively destroy pigment clusters, blood-filled capillary veins, or hair follicle, depending on the intended nature of the treatment.

The treatment

To ensure maximum comfort, your physician should apply cold coupling gel on the treatment area. Depending on your specific condition, you can expect the session to last 15-20 minutes. You may experience a warm sensation as the light is applied to the skin, but the treatment is gentle and should not feel uncomfortable.



Filtered light specifically targets the bacteria in the inflamed sebaceous glands



Intense pulses of light are delivered through the lightguide with integrated cooling to ensure maximum comfort.